

Anti-HMGCS1 Antibody Picoband® APC Conjugated

Catalog Number: A05313-2-APC

About HMGCS1

In molecular biology, hydroxymethylglutaryl-CoA synthase or HMG-CoA synthase EC 2.3.3.10 is an enzyme which catalyzes the reaction in which acetyl-CoA condenses with acetoacetyl-CoA to form 3-hydroxy-3-methylglutaryl-CoA (HMG-CoA). Cytoplasmic 3-hydroxy-3-methylglutaryl-CoA synthase (HMGCS1) catalyzes the first committed step of mevalonate pathway essential for cholesterol biosynthesis. HMGCS1 transcription is regulated by sterol levels. Studies showed that miR-223 reduces cholesterol biosynthesis by inhibiting HMGCS1 and methylsterol monooxygenase 1. In addition, activation of the EGFR family member ERBB4 induces the expression of SREBP-regulated genes (including HMGCS1) involved in cholesterol biosynthesis, suggesting a role of this enzyme in the metabolic re-programming in ERBB4-driven cancers.

Overview

Product Name	Anti-HMGCS1 Antibody Picoband® APC Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, IHC, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q01581

Technical Details

Immunogen	E.coli-derived human HMGCS1 recombinant protein (Position: K251-H520).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HMGCS1 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.