

Anti-KIFC1 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A05325-1-Fluoro594

About KIFC1

Kinesin-like protein KIFC1 is a protein that in humans is encoded by the KIFC1 gene. Kinesin superfamily proteins (KIFs) are molecular motors that drive directional, microtubule-dependent intracellular transport of membrane-bound organelles and other macromolecules (e.g. proteins, nucleic acids). The intracellular transport functions of KIFs are fundamentally important for a variety of cellular functions, including mitotic and meiotic division, motility/migration, hormone and neurotransmitter release, and differentiation. Disruptions to KIF-mediated intracellular transport have been linked with a variety of pathologies, ranging from tumorigenesis to defects in higher order brain function such as learning and memory. KIFC1/HSET is a minus-end directed KIF involved in the processing and movement of early endocytic vesicles, as well as microtubule crosslinking and spindle assembly.

Overview

| | |
|----------------------|--|
| Product Name | Anti-KIFC1 Antibody Picoband® Fluoro594 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q9BW19 |

Technical Details

| | |
|------------------|---|
| Immunogen | E.coli-derived human KIFC1 recombinant protein (Position: M1-K673). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |

Suggested Dilutions

Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-KIFC1 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.