

Anti-KMT5A Antibody

Catalog Number: A05349-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-KMT5A Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-KMT5A Antibody catalog # A05349-1. Tested in WB,Flow Cytometry applications. This antibody reacts with Human,Mouse.
Application	Flow Cytometry, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NQR1


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500-1:2,000 FC: 1:50-1:100 Affinity purified Bach-2 antibody has been tested for use in ELISA, Immunoprecipitation, and western blot. Specific conditions for reactivity should be optimized by the end user. Expect a band approximately 91-140 kDa in size corresponding to Bach2 protein by western blotting in the appropriate cell lysate or extract.

Anti-KMT5A Antibody (A05349-1) Images



Western blot analysis of KMT5A on different lysates using anti-KMT5A antibody at 1/1,000 dilution. Positive control: Lane 1: Mouse kidney Lane 2: Mouse liver Lane 3: 293T Lane 4: 293



Flow cytometric analysis of 293T cells with KMT5A antibody at 1/100 dilution (red) compared with an unlabelled control (cells without incubation with primary antibody; black).

Anti-KMT5A Antibody

For Research Use Only. Not for use in diagnostic procedures.