

Anti-GABP- alpha Monoclonal Antibody

Catalog Number: A05350


Overview

Product Name	Anti-GABP- alpha Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-GABP- alpha Monoclonal Antibody catalog # A05350. Tested in WB, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, WB
Clonality	Monoclonal 8C1B1
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q06546


Technical Details

Immunogen	Purified recombinant fragment of human GABP-alpha (aa120-190) expressed in E. Coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000

Anti-GABP- alpha Monoclonal Antibody (A05350) Images



Confocal immunofluorescence analysis of HeLa cells using GABP-alpha Monoclonal Antibody (green). Red: Actin filaments have been labeled using DY-554 phalloidin. Blue: DRAQ5 fluorescent DNA dye.



Western Blot analysis using GABP-alpha Monoclonal Antibody against HeLa (1), A549 (2), MCF-7 (3), NIH/3T3 (4) and SMMC-7721 (5) cell lysate.

1 Publications Citing This Product

1. PubMed ID: 10.1089/rej.2017.1960, Naringenin Ameliorates Behavioral Dysfunction and Neurological Deficits in a d-Galactose-Induced Aging Mouse Model Through Activation of PI3K/Akt/Nrf2 Pathway

Visit bosterbio.com/anti-gabp-alpha-antibody-a05350-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GABP- alpha Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.