

Anti-Syntrophin Alpha 1 SNTA1 Antibody

Catalog Number: A05370

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Syntrophin Alpha 1 SNTA1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Syntrophin Alpha 1 SNTA1 Antibody catalog # A05370. Tested in WB, ICC, IF, IHC, Flow Cytometry applications. This antibody reacts with Human, Mouse.
Application	Flow Cytometry, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13424


Technical Details

Immunogen	Recombinant protein within human Syntrophin alpha 1 aa 1-200.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 ICC: 1:50-1:200 IHC: 1:50-1:200 FC: 1:50-1:200

Anti-Syntrophin Alpha 1 SNTA1 Antibody (A05370) Images



Western blot analysis of Syntrophin alpha 1 on SiHa cells lysates using anti-Syntrophin alpha 1 antibody at 1/1,000 dilution.



ICC staining Syntrophin alpha 1 in LOVO cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

[Anti-Syntrophin Alpha 1 SNTA1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.