

Anti-Asialoglycoprotein receptor 1 ASGR1 Antibody

Catalog Number: A05376

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Asialoglycoprotein receptor 1 ASGR1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Asialoglycoprotein receptor 1 ASGR1 Antibody catalog # A05376. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P07306

Technical Details

Immunogen	Recombinant fusion protein of human ASGR1(NP_001662.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:1000-1:4000 IHC: 1:50-1:200

Anti-Asialoglycoprotein receptor 1 ASGR1 Antibody (A05376) Images



Western blot analysis of extracts of various cell lines, using ASGR1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.

Anti-Asialoglycoprotein receptor 1 ASGR1 Antibody

For Research Use Only. Not for use in diagnostic procedures.