

Anti-TFIIH p62 GTF2H1 Monoclonal Antibody

Catalog Number: A05408-1


Overview

Product Name	Anti-TFIIH p62 GTF2H1 Monoclonal Antibody
Reactive Species	Bovine, Dog, Human, Mouse, Pig, Rat
Description	Boster Bio Anti-TFIIH p62 GTF2H1 Monoclonal Antibody catalog # A05408-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat, Bovine, Dog, Pig.
Application	WB
Clonality	Monoclonal 8C1B1
Formulation	Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4, 150 mM NaCl) with 0.2% sodium azide, 50% glycerol.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	P32780

Technical Details

Immunogen	Purified recombinant human TFIIH p62 (N-terminus) protein fragments expressed in E.coli.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Protein A purified
Suggested Dilutions	WB 1:1000-1:2000

Anti-TFIIH p62 GTF2H1 Monoclonal Antibody (A05408-1) Images



Western Blot analysis using TFIIH p62 Monoclonal Antibody against LNCAP, K562 cell lysate.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TFIIH p62 GTF2H1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.