

Anti-FGFR-5 FGFR1 Antibody

Catalog Number: A05416


Overview

Product Name	Anti-FGFR-5 FGFR1 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-FGFR-5 FGFR1 Antibody catalog # A05416. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8N441

Technical Details

Immunogen	Synthesized peptide derived from FGFR-5 . at AA range: 130-210
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:40000

Anti-FGFR-5 FGFR-5 Polyclonal Antibody (A05416) Images



Western Blot analysis of 3T3 cells using FGFR-5 Polyclonal Antibody

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FGFR-5 FGFR-5 Polyclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.