

Anti-Endophilin I (P84) SH3GL2 Antibody

Catalog Number: A05430-1


Overview

Product Name	Anti-Endophilin I (P84) SH3GL2 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-Endophilin I (P84) SH3GL2 Antibody catalog # A05430-1. Tested in WB applications. This antibody reacts with Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99962

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 50-99 of Human Endophilin I.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Endophilin I (P84) SH3GL2 Antibody (A05430-1) Images



Western blot (WB) analysis of Endophilin I (P84) pAb at 1:500 dilution Lane1:The Brain tissue lysate of Mouse(10ug) Lane2:The Brain tissue lysate of Rat(10ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Endophilin I (P84) SH3GL2 Antibody

For Research Use Only. Not for use in diagnostic procedures.