

## Anti-HtrA3 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A05478-Fluoro488

### About HTRA3

Human HtrA3 protease, which induces mitochondria-mediated apoptosis, can be a tumor suppressor and a potential therapeutic target in the treatment of cancer. It may also have a role in ovarian development, granulosa cell differentiation and luteinization. The long isoform, HTRA3L, contains 453 amino acids and has a predicted molecular mass of 49 kD. It contains an N-terminal signal peptide, followed by an insulin/IGF (see 147440)-binding domain, a Kazal-type S protease inhibitor domain, a trypsin protease domain, and a PDZ domain. The short isoform, HTRA3S, contains 357 amino acids and has a predicted molecular mass of 38 kD.

### Overview

Product Name	Anti-HtrA3 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P83110

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human HtrA3, different from the related mouse sequence by four amino acids, and from the related rat sequence by three amino acids.
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

---

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HtrA3 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.