

## Anti-CDC42BPA Antibody

Catalog Number: A05509

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




### Overview

|                      |                                                                                                                                            |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-CDC42BPA Antibody                                                                                                                     |
| Reactive Species     | Human, Mouse                                                                                                                               |
| Description          | Boster Bio Anti-CDC42BPA Antibody catalog # A05509. Tested in WB applications. This antibody reacts with Human, Mouse.                     |
| Application          | WB                                                                                                                                         |
| Clonality            | Polyclonal                                                                                                                                 |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2                                                                     |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit                                                                                                                                     |
| Uniprot ID           | Q5VT25                                                                                                                                     |

### Technical Details

|                     |                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen           | A synthetic peptide corresponding to residues in Human CDC42BPA.                                                                                            |
| Isotype             | IgG                                                                                                                                                         |
| Form                | Liquid                                                                                                                                                      |
| Concentration       | 1 mg/ml                                                                                                                                                     |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000                                                                                                                                            |

## Anti-CDC42BPA Antibody (A05509) Images



Western blot (WB) analysis of CDC42BPA pAb at 1:500 dilution Lane1:Panc1 whole cell lysate(40ug) Lane2: U-87MG whole cell lysate(40ug)

### [Anti-CDC42BPA Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.