

Anti-PBX3 Antibody

Catalog Number: A05521-1

About PBX3

Predicted to enable DNA-binding transcription factor activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in animal organ development; neuron development; and regulation of transcription by RNA polymerase II. Predicted to act upstream of or within several processes, including adult locomotory behavior; dorsal spinal cord development; and regulation of respiratory gaseous exchange by nervous system process. Predicted to be located in nucleus. Predicted to be part of chromatin.

Overview

Product Name	Anti-PBX3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PBX3 Antibody catalog # A05521-1. Tested in WB, IHC, ICC, IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	P40426

Technical Details

Immunogen	E.coli-derived human PBX3 recombinant protein (Position: P104-N455).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 Immunoprecipitation, 1:50

ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PBX3 Antibody

For Research Use Only. Not for use in diagnostic procedures.