

Anti-GDF-6 Antibody

Catalog Number: A05536

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-GDF-6 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-GDF-6 Antibody catalog # A05536. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6KF10


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human GDF6. AA range:311-360
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified


Suggested Dilutions

WB 1:500-1:2000
ELISA 1:10000


Anti-GDF-6 Antibody (A05536) Images



Western Blot (WB) analysis of Rat Kidney cells using GDF-6 Polyclonal antibody.



Western Blot (WB) analysis of 293T using GDF-6 Polyclonal antibody diluted at 1:500.



Western Blot (WB) analysis of Rat Kidney cells using GDF-6 Polyclonal antibody.

Anti-GDF-6 Antibody

For Research Use Only. Not for use in diagnostic procedures.