

## Anti-Zic1/2/3 (D348) Antibody

Catalog Number: A05537-1

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview


|                      |                                                                                                                                            |
|----------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-Zic1/2/3 (D348) Antibody                                                                                                              |
| Reactive Species     | Human, Mouse, Rat                                                                                                                          |
| Description          | Boster Bio Anti-Zic1/2/3 (D348) Antibody catalog # A05537-1. Tested in IHC,IF applications. This antibody reacts with Human,Mouse,Rat.     |
| Application          | IF, IHC                                                                                                                                    |
| Clonality            | Polyclonal                                                                                                                                 |
| Formulation          | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2                                                                     |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host                 | Rabbit                                                                                                                                     |
| Uniprot ID           | Q15915/O95409/O60481                                                                                                                       |

### Technical Details

|                     |                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunogen           | Synthetic peptide, corresponding to amino acids 318-366 of Human Zic1.                                                                                      |
| Form                | Liquid                                                                                                                                                      |
| Concentration       | 1 mg/ml                                                                                                                                                     |
| Purification        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | IHC: 1:50-1:200<br>IF: 1:50-1:200                                                                                                                           |

## Anti-Zic1/2/3 (D348) Antibody (A05537-1) Images

---



Immunohistochemistry (IHC) analyzes of Zic1/2/3 (D348) pAb in paraffin-embedded human brain tissue.

[Anti-Zic1/2/3 \(D348\) Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.