

## Anti-MORF4L1 Antibody Picoband®

Catalog Number: A05547-2

### About MORF4L1

Predicted to enable chromatin binding activity. Involved in double-strand break repair via homologous recombination; positive regulation of double-strand break repair via homologous recombination; and regulation of cell cycle. Located in nuclear speck. Part of NuA4 histone acetyltransferase complex and Sin3-type complex.

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Overview


|                      |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-MORF4L1 Antibody Picoband®                                                                                                                                                                                                                                                                                                                                                                                                               |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Description          | Boster Bio Anti-MORF4L1 Antibody Picoband® catalog # A05547-2. Tested in WB, IHC, ICC, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance. |
| Application          | ELISA, IF, IHC, ICC, WB                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Formulation          | Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na <sub>2</sub> HPO <sub>4</sub> .                                                                                                                                                                                                                                                                                                                                                     |
| Storage Instructions | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.                                                                                                                                                                                                                                             |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Uniprot ID           | Q9UBU8                                                                                                                                                                                                                                                                                                                                                                                                                                        |

### Technical Details


|               |                                                                           |
|---------------|---------------------------------------------------------------------------|
| Immunogen     | E.coli-derived human MORF4L1 recombinant protein (Position: M1-E125).     |
| Form          | Lyophilized                                                               |
| Concentration | Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml. |

|                     |                                                                                                                                                                                                             |
|---------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Purification        | Immunogen affinity purified.                                                                                                                                                                                |
| Suggested Dilutions | Western blot, 0.25-0.5 ug/ml, Human, Mouse, Rat<br>Immunohistochemistry(Paraffin-embedded Section), 2-5 ug/ml, Human, Rat<br>Immunocytochemistry/Immunofluorescence, 5 ug/ml, Human<br>ELISA, 0.1-0.5 ug/ml |


## Anti-MORF4L1 Antibody Picoband® (A05547-2) Images



Western blot analysis of MORF4L1 using anti-MORF4L1 antibody (A05547-2). Electrophoresis was performed on a 10% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human A549 whole cell lysates, Lane 2: human MCF-7 whole cell lysates, Lane 3: human Hela whole cell lysates, Lane 4: human HUH-7 whole cell lysates, Lane 5: rat brain tissue lysates, Lane 6: rat liver tissue lysates, Lane 7: mouse brain tissue lysates, Lane 8: mouse liver tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-MORF4L1 antigen affinity purified polyclonal antibody (A05547-2) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for MORF4L1 at approximately 35-38 kDa. The expected band size for MORF4L1 is at 41 kDa.




IHC analysis of MORF4L1 using anti-MORF4L1 antibody (A05547-2). MORF4L1 was detected in a paraffin-embedded section of human skin cancer tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 2 ug/ml rabbit anti-MORF4L1 Antibody (A05547-2) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.



IHC analysis of MORF4L1 using anti-MORF4L1 antibody (A05547-2). MORF4L1 was detected in a paraffin-embedded section of rat testis tissue. Heat mediated antigen retrieval was performed in EDTA buffer (pH 8.0, epitope retrieval solution). The tissue section was blocked with 10% goat serum. The tissue section was then incubated with 2 ug/ml rabbit anti-MORF4L1 Antibody (A05547-2) overnight at 4°C. Peroxidase Conjugated Goat Anti-rabbit IgG was used as secondary antibody and incubated for 30 minutes at 37°C. The tissue section was developed using HRP Conjugated Rabbit IgG Super Vision Assay Kit (Catalog # SV0002) with DAB as the chromogen.

IF analysis of MORF4L1 using anti-MORF4L1 antibody (A05547-2) and anti-Alpha Tubulin antibody (M03989-3).



MORF4L1 was detected in an immunocytochemical section of HeLa cells. Enzyme antigen retrieval was performed using IHC enzyme antigen retrieval reagent (AR0022) for 15 mins. The cells were blocked with 10% goat serum. And then incubated with 5 ug/mL rabbit anti-MORF4L1 Antibody (A05547-2) and mouse anti-Alpha Tubulin antibody (M03989-3) overnight at 4°C. Cy3 Conjugated Goat Anti-Rabbit IgG (BA1032) and Fluoro488 Conjugated Goat Anti-Mouse IgG (BA1126) were used as secondary antibody at 1:500 dilution and incubated for 30 minutes at 37°C. The section was counterstained with DAPI. Visualize using a fluorescence microscope and filter sets appropriate for the label used.

#### [Anti-MORF4L1 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.