

## Anti-JunD (V249) Antibody

Catalog Number: A05609


### Overview

Product Name	Anti-JunD (V249) Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-JunD (V249) Antibody catalog # A05609. Tested in WB,IHC,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P17535

### Technical Details


Immunogen	Synthetic peptide, corresponding to amino acids 220-270 of Human JunD.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IP: 1:50-1:200

## Anti-JunD (V249) Antibody (A05609) Images



A549 whole cell lysate  
JunD (V249) pAb at 1:500 dilution

Western blot (WB) analysis of JunD (V249) polyclonal antibody in extracts from A549 cells.



Immunohistochemistry (IHC) analyzes of JunD (V249) pAb in paraffin-embedded human breast carcinoma tissue at 1:100.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-JunD (V249) Antibody

For Research Use Only. Not for use in diagnostic procedures.