

Anti-Plakophilin 4 PKP4 Antibody

Catalog Number: A05639

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Plakophilin 4 PKP4 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Plakophilin 4 PKP4 Antibody catalog # A05639. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q99569


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human PKP4. AA range:185-234
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified

Suggested Dilutions

WB 1:500-1:2000
IHC 1:100-1:300
ELISA 1:5000
IF 1:50-200

Anti-Plakophilin 4 PKP4 Antibody (A05639) Images



Western Blot (WB) analysis of Mouse Heart lysis using Plakophilin 4 antibody.

Anti-Plakophilin 4 PKP4 Antibody

For Research Use Only. Not for use in diagnostic procedures.