

Anti-Frizzled-2 FZD2 Antibody

Catalog Number: A05675-1

Overview

Product Name	Anti-Frizzled-2 FZD2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Frizzled-2 FZD2 Antibody catalog # A05675-1. Tested in ELISA, IF, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14332

Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human FZD2. AA range:201-250
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:40000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Frizzled-2 FZD2 Antibody