

Anti-MLL4 KMT2B Antibody

Catalog Number: A05680-1

Overview

Product Name	Anti-MLL4 KMT2B Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-MLL4 KMT2B Antibody catalog # A05680-1. Tested in IHC, IF applications. This antibody reacts with Human, Mouse.
Application	IF, IHC
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UMN6

Technical Details

Immunogen	Synthesized peptide derived from human protein at AA range: 1060-1140
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	MLL4 antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	IHC 1:50-300 IF 1:50-200

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MLL4 KMT2B Antibody

For Research Use Only. Not for use in diagnostic procedures.