

## Anti-PIPK I gamma PIP5K1C Antibody

Catalog Number: A05687


### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-PIPK I gamma PIP5K1C Antibody   |
| Reactive Species     | Human, Mouse, Rat  |
| Description          | Boster Bio Anti-PIPK I gamma PIP5K1C Antibody catalog # A05687. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.   |
| Application          | ELISA, WB  |
| Clonality            | Polyclonal   |
| Formulation          | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide.<br>*This antibody is supplied in a stabilized formulation.<br>Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.   |
| Host                 | Rabbit   |
| Uniprot ID           | O60331   |


### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | The antiserum was produced against synthesized peptide derived from human PIP5K1C. AA range:305-354 |
| Cross Reactivity    | No cross reactivity with other proteins.  |
| Isotype             | IgG   |
| Form                | Liquid  |
| Concentration       | 1 mg/ml   |
| Purification        | Immunogen affinity purified   |
| Suggested Dilutions | WB 1:500-1:2000<br>ELISA 1:20000  |


## Anti-PIP5K1C Antibody (A05687) Images



Western Blot analysis of A549 cells using PIPK I gamma Polyclonal Antibody



Western blot analysis of lysates from A549 cells, using PIP5K1C Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HepG2 cells using PIP5K1C antibody.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PIP5K1C Antibody

For Research Use Only. Not for use in diagnostic procedures.