

Anti-PIP5KI Gamma PIP5K1C Antibody

Catalog Number: A05687-1

Overview

Product Name	Anti-PIP5KI Gamma PIP5K1C Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-PIP5KI Gamma PIP5K1C Antibody catalog # A05687-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O60331

Technical Details

Form	Liquid
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	<p>Dilute the sample so that the expected range of concentrations fall within the detection range of this kit.</p> <p>If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples.</p> <p>Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB: 1:500~1:1000</p> <p>For protocols, please visit https://www.bosterbio.com/protocol-and-troubleshooting/</p>

Anti-PIP5KI Gamma PIP5K1C Antibody (A05687-1) Images




Figure 1. Western blotting validation for Anti-PIP5KI Gamma PIP5K1C Antibody A05687-1

Western blot (WB) analysis of PIP5K1 γ polyclonal antibody at 1:500 dilution

Lane1:HeLa cell lysate

Lane2:Raw264.7 cell lysate

Lane3:PC12 cell lysate

Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PIP5KI Gamma PIP5K1C Antibody