

## Anti-eIF5A2 (Y127) Antibody

Catalog Number: A05688Y127


### Overview

Product Name	Anti-eIF5A2 (Y127) Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-eIF5A2 (Y127) Antibody catalog # A05688Y127. Tested in WB,IHC applications. This antibody reacts with Human,Mouse.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9GZV4

### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 96-140 of Human eIF5A2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

## Anti-eIF5A2 (Y127) Antibody (A05688Y127) Images



Western blot (WB) analysis of eIF5A2 (Y127) pAb at 1:500 dilution Lane1:AML-12 whole cell lysate(40ug) Lane2:C6 whole cell lysate(40ug) Lane3:HEK293T whole cell lysate(40ug) Lane4:SGC7901 whole cell lysate(40ug) Lane5:HepG2 whole cell lysate(40ug) Lane6:A549 whole cell lysate(40ug) Lane7:A2780 whole cell lysate(40ug)

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-eIF5A2 (Y127) Antibody

For Research Use Only. Not for use in diagnostic procedures.