

Anti-Ras-related protein Rap-2b RAP2B Antibody

Catalog Number: A05702


Overview

Product Name	Anti-Ras-related protein Rap-2b RAP2B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Ras-related protein Rap-2b RAP2B Antibody catalog # A05702. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P61225

Technical Details

Immunogen	Synthetic peptide, corresponding Human RAP2B.
Form	Liquid
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-Ras-related protein Rap-2b RAP2B Antibody (A05702) Images



Western blot (WB) analysis of RAP2B polyclonal antibody at 1:500 dilution Lane1:K562 whole cell lysate(40ug) Lane2:Jurkat whole cell lysate(40ug) Lane3:THP-1 whole cell lysate(40ug) Lane4:BV2 whole cell lysate(40ug) Lane5:PMVEC whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ras-related protein Rap-2b RAP2B Antibody

For Research Use Only. Not for use in diagnostic procedures.