

## Anti-TID1/DNAJA3 Antibody Picoband® FITC Conjugated

Catalog Number: A05761-1-FITC

### About DNAJA3

This gene encodes a member of the DNAJ/Hsp40 protein family. DNAJ/Hsp40 proteins stimulate the ATPase activity of Hsp70 chaperones and play critical roles in protein folding, degradation, and multimeric complex assembly. The encoded protein is localized to mitochondria and mediates several cellular processes including proliferation, survival and apoptotic signal transduction. The encoded protein also plays a critical role in tumor suppression through interactions with oncogenic proteins including ErbB2 and the p53 tumor suppressor protein. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-TID1/DNAJA3 Antibody Picoband® FITC Conjugated  |
| Reactive Species     | Human  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | Q96EY1   |

### Technical Details

|                     |  |
|---------------------|--|
| Immunogen           | E.coli-derived human TID1/DNAJA3 recombinant protein (Position: K73-G436). |
| Form                | Liquid   |
| Concentration       | 0.5 mg/mL  |
| Purification        | Immunogen affinity purified.   |
| Conjugate           | FITC<br>Excitation Wavelength: 495 nm<br>Emission Wavelength: 525 nm       |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.       |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-TID1/DNAJA3 Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.