

## Anti-SASS6 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A05762-1-Fluoro550

### About SASS6

Spindle assembly abnormal protein 6 homolog (SAS-6) is a protein that in humans is encoded by the SASS6 gene. The protein encoded by this gene is a central component of centrioles and is necessary for their duplication and function. Centrioles adopt a cartwheel-shaped structure, with the encoded protein forming the hub and spokes inside a microtubule cylinder. Defects in this gene are a cause of autosomal recessive primary microcephaly.

### Overview

Product Name	Anti-SASS6 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6UVJ0

### Technical Details

Immunogen	E.coli-derived human SASS6 recombinant protein (Position: E101-S657). Human SASS6 shares 83.7% amino acid (aa) sequence identity with mouse SASS6.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SASS6 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.