

Anti-Pseudokinase FAM20A FAM20A Antibody

Catalog Number: A05777

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Pseudokinase FAM20A FAM20A Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-Pseudokinase FAM20A FAM20A Antibody catalog # A05777. Tested in WB applications. This antibody reacts with Human, Mouse.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q96MK3


Technical Details

Immunogen	KLH conjugated synthetic peptide derived from human FAM20A:101-200/541
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-Pseudokinase FAM20A FAM20A Antibody (A05777) Images



Sample: Lung(Mouse) Lysate at 40 ug Primary: Anti-FA20A
at 1/300 dilution



Sample: MDA-MB-231 (human)Cell Lysate at 40 ug Primary:
Anti-FA20A at 1/300 dilution

[Anti-Pseudokinase FAM20A FAM20A Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.