

Anti-CDKN2AIP/CARF Antibody Picoband®

Catalog Number: A05786-1

About Carf

Enables DNA-binding transcription activator activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Involved in positive regulation of transcription by RNA polymerase II. Acts upstream of or within regulation of DNA-binding transcription factor activity. Located in granular component. Is expressed in several structures, including central nervous system; eye; genitourinary system; gut; and skeletal muscle tissue. Orthologous to human CARF (calcium responsive transcription factor).

Overview

Product Name	Anti-CDKN2AIP/CARF Antibody Picoband®
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-CDKN2AIP-CARF Antibody Picoband® catalog # A05786-1. Tested in WB, Flow Cytometry, ELISA applications. This antibody reacts with Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, Flow Cytometry, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8VHI4


Technical Details

Immunogen	E.coli-derived mouse CDKN2AIP/CARF recombinant protein (Position: R183-T689). Human CDKN2AIP/CARF shares 78.3% and 88.1% amino acid (aa) sequence identity with human and rat CDKN2AIP/CARF, respectively.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.


Suggested Dilutions

Western blot, 0.25-0.5 ug/ml, Mouse, Rat
Flow Cytometry (Fixed), 1-3 ug/ 1×10^6 cells, Mouse
ELISA, 0.1-0.5 ug/ml

Anti-CDKN2AIP/CARF Antibody Picoband® (A05786-1) Images



Western blot analysis of CDKN2AIP/CARF using anti-CDKN2AIP/CARF antibody (A05786-1). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: rat brain tissue lysates, Lane 2: rat lung tissue lysates, Lane 3: rat kidney tissue lysates, Lane 4: mouse brain tissue lysates, Lane 5: mouse lung tissue lysates, Lane 6: mouse kidney tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-CDKN2AIP/CARF antigen affinity purified polyclonal antibody (A05786-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for CDKN2AIP/CARF at approximately 75 kDa. The expected band size for CDKN2AIP/CARF is at 61 kDa.



Flow Cytometry analysis of NIH/3T3 cells using anti-CDKN2AIP/CARF antibody (A05786-1). Overlay histogram showing NIH/3T3 cells stained with A05786-1 (Blue line). To facilitate intracellular staining, cells were fixed with 4% paraformaldehyde and permeabilized with permeabilization buffer. The cells were blocked with 10% normal goat serum. And then incubated with rabbit anti-CDKN2AIP/CARF Antibody (A05786-1, 1 ug/1x10⁶ cells) for 30 min at 20°C. DyLight® 488 conjugated goat anti-rabbit IgG (BA1127, 5-10 ug/1x10⁶ cells) was used as secondary antibody for 30 minutes at 20°C. Isotype control antibody (Green line) was rabbit IgG (1 ug/1x10⁶) used under the same conditions. Unlabelled sample without incubation with primary antibody and secondary antibody (Red line) was used as a blank control.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CDKN2AIP/CARF Antibody

For Research Use Only. Not for use in diagnostic procedures.