

Anti-Tau-tubulin kinase 2 TTBK2 Antibody

Catalog Number: A05821-1


Overview

Product Name	Anti-Tau-tubulin kinase 2 TTBK2 Antibody
Reactive Species	Human, Rat
Description	Boster Bio Anti-Tau-tubulin kinase 2 TTBK2 Antibody catalog # A05821-1. Tested in WB,IHC applications. This antibody reacts with Human,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q6IQ55


Technical Details

Immunogen	Recombinant protein of human TTBK2.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC 1:50-1:200

Anti-Tau-tubulin kinase 2 TTBK2 Antibody (A05821-1) Images



Western blot analysis of extracts of rat brain, using TTBK2 antibody.



Immunohistochemistry of paraffin-embedded human gastric using TTBK2 antibody at dilution of 1:200 (40x lens).

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Tau-tubulin kinase 2 TTBK2 Antibody

For Research Use Only. Not for use in diagnostic procedures.