

Anti-p70 S6 kinase beta RPS6KB2 Antibody

Catalog Number: A05845


Overview

Product Name	Anti-p70 S6 kinase beta RPS6KB2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-p70 S6 kinase beta RPS6KB2 Antibody catalog # A05845. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UBS0


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human p70 S6 Kinase beta. AA range:389-438
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:5000 IF 1:50-200


Anti-p70 S6 kinase beta RPS6KB2 Antibody (A05845) Images




Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).




Western Blot analysis of various cells using p70 S6 kinase beta Polyclonal Antibody



Western Blot analysis of 293 cells using p70 S6 kinase beta Polyclonal Antibody



Western blot analysis of lysates from 293 cells, using p70 S6 Kinase beta Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HUVEC cells using p70 S6 Kinase beta antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p70 S6 kinase beta RPS6KB2 Antibody

For Research Use Only. Not for use in diagnostic procedures.