

Anti-p70 S6K Beta RPS6KB2 Antibody

Catalog Number: A05845-1


Overview

| | |
|----------------------|--|
| Product Name | Anti-p70 S6K Beta RPS6KB2 Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-p70 S6K Beta RPS6KB2 Antibody catalog # A05845-1. Tested in WB,ICC,IHC,Flow Cytometry applications. This antibody reacts with Human,mouse,rat. |
| Application | Flow Cytometry, IHC, ICC, WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q9UBS0 |


Technical Details

| | |
|---------------------|--|
| Immunogen | peptide |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Peptide affinity purified |
| Suggested Dilutions | WB: 1:500 ICC: 1:200 IHC: 1:200 FC: 1:100-1:200 |

Anti-p70 S6K Beta RPS6KB2 Antibody (A05845-1) Images



Western blot analysis of p70 S6 Kinase beta on different cell lysates using anti-p70 S6 Kinase beta antibody at 1/500 dilution. Positive control: Lane 1: PC12 Lane 2: NIH/3T3 Lane 3: Jurkat Lane 4: K562 Lane 5: Hela



ICC staining p70 S6 Kinase beta in Lovo cells (green). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-p70 S6K Beta RPS6KB2 Antibody

For Research Use Only. Not for use in diagnostic procedures.