

## Anti-STRN Antibody

Catalog Number: A05855-2

### About STRN

Enables armadillo repeat domain binding activity; nuclear estrogen receptor binding activity; and protein phosphatase 2A binding activity. Involved in Wnt signaling pathway and negative regulation of cell population proliferation. Located in bicellular tight junction. Part of FAR/SIN/STRIPAK complex.

### Overview

Product Name	Anti-STRN Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-STRN Antibody catalog # A05855-2. Tested in WB, IHC, ICC/IF, IP, Flow Cytometry, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, Flow Cytometry, IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol *This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	O43815

### Technical Details

Immunogen	E.coli-derived human STRN recombinant protein (Position: 303-400).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400 Immunocytochemistry/Immunofluorescence, 1:50-400 Immunoprecipitation, 1:250-300 Flow Cytometry (Fixed), 1:50-200

ELISA, 1:100-1000

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-STRN Antibody

For Research Use Only. Not for use in diagnostic procedures.