

Anti-NDUFS3 Antibody

Catalog Number: A05867-1


Overview

Product Name	Anti-NDUFS3 Antibody
Reactive Species	Human
Description	Boster Bio Anti-NDUFS3 Antibody catalog # A05867-1. Tested in WB,IHC,ICC,IF,IP applications. This antibody reacts with Human.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O75489


Technical Details

Immunogen	Recombinant protein
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	ProA affinity purified
Suggested Dilutions	WB: 1:500 IHC: 1:50 ICC: 1:50 IP: 1:10-1:50

Anti-NDUFS3 Antibody (A05867-1) Images



Western blot analysis of NDUFS3 on HepG2 cell using anti-NDUFS3 antibody at 1/500 dilution.



ICC staining NDUFS3 in MCF-7 cells (green). The nuclear counter stain is DAPI (blue). Cells were fixed in paraformaldehyde, permeabilised with 0.25% Triton X100/PBS.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NDUFS3 Antibody

For Research Use Only. Not for use in diagnostic procedures.