

Anti-HSPC150/UBE2T Antibody Picoband® APC Conjugated

Catalog Number: A05874-3-APC

About UBE2T

The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-HSPC150/UBE2T Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NPD8

Technical Details

Immunogen	E.coli-derived human HSPC150/UBE2T recombinant protein (Position: R9-Q145). Human HSPC150/UBE2T shares 90.5% amino acid (aa) sequence identity with mouse HSPC150/UBE2T.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-HSPC150/UBE2T Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.