

Anti-HSPC150/UBE2T Antibody Picoband® FITC Conjugated

Catalog Number: A05874-3-FITC

About UBE2T

The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.

Overview

| | |
|----------------------|--|
| Product Name | Anti-HSPC150/UBE2T Antibody Picoband® FITC Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q9NPD8 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human HSPC150/UBE2T recombinant protein (Position: R9-Q145). Human HSPC150/UBE2T shares 90.5% amino acid (aa) sequence identity with mouse HSPC150/UBE2T. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | FITC Excitation Wavelength: 495 nm Emission Wavelength: 525 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-HSPC150/UBE2T Antibody - FITC

For Research Use Only. Not for use in diagnostic procedures.