

Anti-HSPC150/UBE2T Antibody Picoband® HRP Conjugated

Catalog Number: A05874-3-HRP

About UBE2T

The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.

Overview

Product Name	Anti-HSPC150/UBE2T Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NPD8

Technical Details

Immunogen	E.coli-derived human HSPC150/UBE2T recombinant protein (Position: R9-Q145). Human HSPC150/UBE2T shares 90.5% amino acid (aa) sequence identity with mouse HSPC150/UBE2T.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HSPC150/UBE2T Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.