

Anti-NDP52/CALCOCO2 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A05876-1-Fluoro594

About CALCOCO2

This gene encodes a coiled-coil domain-containing protein. The encoded protein functions as a receptor for ubiquitin-coated bacteria and plays an important role in innate immunity by mediating macroautophagy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

Overview

| | |
|----------------------|--|
| Product Name | Anti-NDP52/CALCOCO2 Antibody Picoband® Fluoro594 Conjugated |
| Reactive Species | Human |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q13137 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human NDP52/CALCOCO2 recombinant protein (Position: M1-L446). |
| Cross Reactivity | No cross-reactivity with other proteins. |
| Isotype | Rabbit IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NDP52/CALCOCO2 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.