

Anti-Ras-related protein Rab-31 RAB31 Antibody

Catalog Number: A05883-1


Overview

Product Name	Anti-Ras-related protein Rab-31 RAB31 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Ras-related protein Rab-31 RAB31 Antibody catalog # A05883-1. Tested in WB, ICC applications. This antibody reacts with Human, Mouse, Rat.
Application	ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q13636

Technical Details

Immunogen	Synthetic peptide, corresponding to the 121-166 of Human RAB31.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 ICC: 1:50-1:200

Anti-Ras-related protein Rab-31 RAB31 Antibody (A05883-1) Images



Western blot (WB) analysis of RAB31 polyclonal antibody at 1:500 dilution Lane1:HEK293T whole cell lysate(40ug) Lane2:U-87MG whole cell lysate(40ug) Lane3:C6 whole cell lysate(40ug) Lane4:CT26 whole cell lysate(40ug) Lane5:L02 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Ras-related protein Rab-31 RAB31 Antibody

For Research Use Only. Not for use in diagnostic procedures.