

Anti-RAB31 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A05883-2-Fluoro594

About RAB31

Enables GDP binding activity; GTP binding activity; and GTPase activity. Involved in several processes, including Golgi to plasma membrane protein transport; cellular response to insulin stimulus; and receptor internalization. Located in early endosome; phagocytic vesicle; and trans-Golgi network membrane. Biomarker of severe acute respiratory syndrome.

Overview

| | |
|----------------------|---|
| Product Name | Anti-RAB31 Antibody Picoband® Fluoro594 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q13636 |

Technical Details

| | |
|---------------------|---|
| Immunogen | A synthetic peptide corresponding to a sequence at the C-terminus of human RAB31. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAB31 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.