

Anti-RAB31 Antibody Picoband® HRP Conjugated

Catalog Number: A05883-2-HRP

About RAB31

Enables GDP binding activity; GTP binding activity; and GTPase activity. Involved in several processes, including Golgi to plasma membrane protein transport; cellular response to insulin stimulus; and receptor internalization. Located in early endosome; phagocytic vesicle; and trans-Golgi network membrane. Biomarker of severe acute respiratory syndrome.

Overview

Product Name	Anti-RAB31 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q13636

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the C-terminus of human RAB31.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RAB31 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.