

## Anti-SSX2IP Antibody Picoband® Biotin Conjugated

Catalog Number: A05918-1-Biotin

### About SSX2IP

Afadin- and alpha-actinin-binding protein is a protein that in humans is encoded by the SSX2IP gene. This gene encodes a protein that binds the cancer-testis antigen Synovial Sarcoma X breakpoint 2 protein. The encoded protein may regulate the activity of Synovial Sarcoma X breakpoint 2 protein in malignant cells. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 3.

### Overview

Product Name	Anti-SSX2IP Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y2D8

### Technical Details

Immunogen	E.coli-derived human SSX2IP recombinant protein (Position: K27-P614). Human SSX2IP shares 87.9% and 87.4% amino acid (aa) sequence identity with mouse and rat SSX2IP, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-SSX2IP Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.