

Anti-SSX2IP Antibody Picoband® Cy3 Conjugated

Catalog Number: A05918-1-Cy3

About SSX2IP

Afadin- and alpha-actinin-binding protein is a protein that in humans is encoded by the SSX2IP gene. This gene encodes a protein that binds the cancer-testis antigen Synovial Sarcoma X breakpoint 2 protein. The encoded protein may regulate the activity of Synovial Sarcoma X breakpoint 2 protein in malignant cells. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 3.

Overview

| | |
|----------------------|--|
| Product Name | Anti-SSX2IP Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human, Rat |
| Application | Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs. |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q9Y2D8 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human SSX2IP recombinant protein (Position: K27-P614). Human SSX2IP shares 87.9% and 87.4% amino acid (aa) sequence identity with mouse and rat SSX2IP, respectively. |
| Cross Reactivity | No cross reactivity with other proteins. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Optimal dilutions should be determined by end users. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SSX2IP Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.