

## Anti-SSX2IP Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A05918-1-Fluoro594

### About SSX2IP

Afadin- and alpha-actinin-binding protein is a protein that in humans is encoded by the SSX2IP gene. This gene encodes a protein that binds the cancer-testis antigen Synovial Sarcoma X breakpoint 2 protein. The encoded protein may regulate the activity of Synovial Sarcoma X breakpoint 2 protein in malignant cells. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 3.

### Overview

Product Name	Anti-SSX2IP Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9Y2D8

### Technical Details

Immunogen	E.coli-derived human SSX2IP recombinant protein (Position: K27-P614). Human SSX2IP shares 87.9% and 87.4% amino acid (aa) sequence identity with mouse and rat SSX2IP, respectively.
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SSX2IP Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.