

## Anti-PELI1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A05921-1-Fluoro488

### About PELI1

Enables ubiquitin-ubiquitin ligase activity. Involved in several processes, including negative regulation of TORC1 signaling; negative regulation of necroptotic process; and protein polyubiquitination. Is active in site of double-strand break. Implicated in colorectal cancer.

### Overview

Product Name	Anti-PELI1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (Flow Cytometry, IP, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q96FA3

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PELI1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-PELI1 Antibody - Fluoro488**

**For Research Use Only. Not for use in diagnostic procedures.**