

Anti-PELI1 Antibody Picoband® HRP Conjugated

Catalog Number: A05921-1-HRP

About PELI1

Enables ubiquitin-ubiquitin ligase activity. Involved in several processes, including negative regulation of TORC1 signaling; negative regulation of necroptotic process; and protein polyubiquitination. Is active in site of double-strand break. Implicated in colorectal cancer.

Overview

Product Name	Anti-PELI1 Antibody Picoband® HRP Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q96FA3

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PELI1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	HRP
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



For Research Use Only. Not for use in diagnostic procedures.