

## Anti-NPT1 SLC17A1 Monoclonal Antibody

Catalog Number: A05923


### Overview

Product Name	Anti-NPT1 SLC17A1 Monoclonal Antibody
Reactive Species	Human
Description	Boster Bio Anti-NPT1 SLC17A1 Monoclonal Antibody catalog # A05923. Tested in WB, ELISA applications. This antibody reacts with Human.
Application	ELISA, WB
Clonality	Monoclonal 4B4D1
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q14916

### Technical Details

Immunogen	Purified recombinant fragment of NPT1 expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:10000

## Anti-NPT1 SLC17A1 Monoclonal Antibody (A05923) Images



Western Blot (WB) analysis using NPT1 Monoclonal antibody against truncated NPT1 recombinant protein.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-NPT1 SLC17A1 Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.