

## Anti-Son of sevenless homolog 2 SOS2 Antibody

Catalog Number: A05957-1

### Overview


Product Name	Anti-Son of sevenless homolog 2 SOS2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Son of sevenless homolog 2 SOS2 Antibody catalog # A05957-1. Tested in ICC, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q07890

### Technical Details

Form	Liquid
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows: Boster Bio's internal QC testing used WB: 1:500~1:1000 ICC: 1:50~1:200  For protocols, please visit <a href="https://www.bosterbio.com/protocol-and-troubleshooting/">https://www.bosterbio.com/protocol-and-troubleshooting/</a>

### Anti-Son of sevenless homolog 2 SOS2 Antibody (A05957-1) Images

Figure 1. Western blotting validation for Anti-Son of sevenless homolog 2 SOS2 Antibody A05957-1



Western blot (WB) analysis of SOS2 polyclonal antibody at 1:500 dilution  
Lane1:Hela cell lysate  
Lane2:sp2/0 cell lysate  
Lane3:H9C2 cell lysate  
Electrophoresis was performed on a SDS-PAGE gel. To determine SDS-PAGE gel concentration

## Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Son of sevenless homolog 2 SOS2 Antibody