

Anti-Collagen alpha-3(IX) chain COL9A3 Antibody

Catalog Number: A05963-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Collagen alpha-3(IX) chain COL9A3 Antibody
Reactive Species	Mouse
Description	Boster Bio Anti-Collagen alpha-3 (IX) chain COL9A3 Antibody catalog # A05963-1. Tested in IHC applications. This antibody reacts with Mouse.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14050

Technical Details

Immunogen	Recombinant fusion protein of human COL9A3(NP_001844.3).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:100-1:200

Anti-Collagen alpha-3(IX) chain COL9A3 Antibody (A05963-1) Images



Immunohistochemistry of paraffin-embedded mouse heart using COL9A3 antibody at dilution of 1:100 .Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

[Anti-Collagen alpha-3\(IX\) chain COL9A3 Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.