

Anti-STAM-binding protein STAMBP Antibody

Catalog Number: A05964

About STAMBP

Promotes the exchange of Ras-bound GDP by GTP By similarity.

Mural R.J., Submitted (SEP-2005) to the EMBL/GenBank/DDBJ databases.

The MGC Project Team, Genome Res. 14:2121-2127(2004).

Overview

Product Name	Anti-STAM-binding protein STAMBP Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-STAM-binding protein STAMBP Antibody catalog # A05964. Tested in WB,IHC,ICC/IF,IP applications. This antibody reacts with Human,Mouse,Rat.
Application	IP, IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O95630

Technical Details

Immunogen	Synthesized peptide derived from human NKp44.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:
WB: 1:500-1:2000
IHC: 1:50-1:200
ICC/IF: 1:50-1:100
IP: 1:20-1:50

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-STAM-binding protein STAMBP Antibody