

Anti-INP5E INPP5E Antibody

Catalog Number: A05973

About INPP5E

Promotes the exchange of Ras-bound GDP by GTP By similarity.

Mural R.J., Submitted (SEP-2005) to the EMBL/GenBank/DDBJ databases.
The MGC Project Team, Genome Res. 14:2121-2127(2004).

Overview

Product Name	Anti-INP5E INPP5E Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-INP5E INPP5E Antibody catalog # A05973. Tested in ELISA, WB applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid form in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NRR6

Technical Details

Immunogen	Synthesized peptide derived from human INP5E protein.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	INP5E antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Suggested Dilutions	Dilute the sample so that the expected range of concentrations fall within the detection range of this kit. If the expected range of concentration is unknown, a pilot test should be conducted to decide the optimal dilution ratio for your samples. Some PubMed article(s) citing the expression level of this target are as follows:

Boster Bio's internal QC testing used:
WB 1:500-2000 ELISA 1:5000-20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-INP5E INPP5E Antibody