

Anti-HAPLN1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A05980-3-Fluoro488

About HAPLN1

Predicted to enable hyaluronic acid binding activity. Predicted to be an extracellular matrix structural constituent conferring compression resistance. Predicted to be involved in nervous system development; positive regulation of neuroblast proliferation; and skeletal system development. Located in collagen-containing extracellular matrix.

Overview

Product Name	Anti-HAPLN1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P10915

Technical Details

Immunogen	E.coli-derived human HAPLN1 recombinant protein (Position: D16-N354). Human HAPLN1 shares 96.5% amino acid (aa) sequence identity with both mouse and rat HAPLN1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-HAPLN1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.